GORDON
ELLIS
& CO

PRIMA[®] GRAB BAR

ANGLED & STRAIGHT

INSTRUCTIONS FOR USE **en**

cs

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da

BRUGSANVISNING
7

de

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es

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fi

KÄYTTÖOHJE
14

fr

MODE D'EMPLOI
16

hu

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nl

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no

BRUKSANVISNING
22

pl

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

sk

NÁVOD NA POUŽITIE
26

sv

BRUKSANVISNING
28

REF

6519, 6520, 6521, 6522, 6523, 6524, 6538, 6540, 6541, 6542, 6553, 6554, 6555,
6556, 6557, 6558

EU REP

Advena Ltd.

Tower Business Centre, 2nd Floor, Tower Street, Swatar. MALTA. BKR 4013

Gordon Ellis & Co.

Trent Lane, Castle Donington, Derbys. UK. DE74 2AT

6519-0926-7

UK
CA

MD

ISS DATE: SEPT 2026

GORDONELLISHEALTHCARE.CO.UK

en Accessible format available on request.
bg Достъпен формат е наличен при поискване.
cs Přístupný formát je k dispozici na vyžádání.
da Tilgængeligt format fås efter anmodning.
de Barrierefreies Format auf Anfrage erhältlich.
el Διαθέσιμη προσβάσιμη μορφή κατόπιν αιτήματος.
es Formato accesible disponible a petición.
et Juurdepääsetav vorming on saadaval nõudmisel.
fi Saavutettava versio saatavilla pyynnöstä.
fr Format accessible disponible sur demande.
hr Dostupan format na zahtjev.


hu Hozzáférhető formátum kérésre elérhető.
it Formato accessibile disponibile su richiesta.
lt Prieinamas formatas pagal pareikalavimą.
lv Pieejams pielāgots formāts pēc pieprasījuma.
nl Toegankelijk formaat op verzoek beschikbaar.
no Tilgjengelig format tilgjengelig på forespørsel.
pl Dostępny format na życzenie.
pt Formato acessível disponível mediante pedido.
ro Format accesibil disponibil la cerere.
sk Dostupný prístupný formát na požiadanie.
sl Na zahtevo je na voljo dostopna oblika.
sv Tillgängligt format finns på begäran.

1


a.

- 1.
- 2.
- 3.

b.




c.



2

12" STRAIGHT	16" STRAIGHT	18" STRAIGHT	24" STRAIGHT	13" ANGLED	16" ANGLED
6521; 6540; 6549; 6555	6522; 6541; 6550; 6556	6523; 6542; 6551; 6557; 6537; 6538	6524; 6552; 6558	6519; 6547; 6553	6520; 6548; 6554

1. 	37 mm					
2. 	177 mm	284 mm	327 mm	477 mm	214 mm	370 mm
3. 	37 x 31 mm					
4. 	L342 x W65 x H68.5 mm	L449 x W65 x H68.5 mm	L492 x W65 x H68.5 mm	L642 x W65 x H68.5 mm	L329 x W65 x H68.5 mm	L457 x W65 x H68.5 mm
5. A 	A = 28.5, B = 13 mm. Ø5.5					


1.


2.


3.


4.

5.

6.

7.

PH 5 - 4

Prior to reading these instructions for use, please check you have the latest version of this document, you can check the latest PDF version online at: <https://www.gordonellishealthcare.co.uk/prima-grab-bars-2>

Here you will also find other product documentation.

Introduction

Welcome to your new Gordon Ellis product. We hope that it proves a useful aid in your daily life. Our UK design team have developed your Prima® Grab Bar specifically with the intention of compensating for an injury or overcoming a temporary or permanent disability. Hopefully making daily living that little bit easier.

These instructions for use contain important information which should be followed to ensure your safety while using this product. Please ensure these instructions are saved for future reference.

Any serious accident or adverse event occurring through use of the device should be reported in a timely manner to your local distributor and the competent authority in your country. The local distributor will forward information to the manufacturer.

Device description

The Prima® Grab Bar range comprises rigid, wall-mounted safety rails designed to provide stable physical support, balance assistance, and leverage for individuals with reduced mobility, physical weakness, or temporary or permanent disabilities.

The range includes both straight and angled configurations across multiple lengths to accommodate varied installation layouts, transfers, and functional user requirements:

Straight Variants: Available in 12" (300 mm), 16" (400 mm), 18" (450 mm), and 24" (600 mm) nominal lengths, engineered for horizontal, or vertical, mounting.

Angled Variants: Available in 13" (330 mm) and 16" (400 mm) nominal lengths, featuring a fixed 130° cranked profile designed to position one section horizontally for downward leverage and pushing up, while the angled section provides support for pulling up and positional transfers.

Each grab bar is constructed around an internal, high-strength aluminium core encased in a moulded plastic exterior. The gripping profile features an ergonomic oval cross-section (37 × 31 mm) integrated with non-slip, ribbed soft-grip inserts that provide secure tactile purchase and remain warm to the touch even in wet bathroom environments.

The assemblies feature concealed wall mounting plates at each terminal flange. Dedicated push-down end caps fully enclose the fixing points, eliminating exposed fastener heads, sharp edges, and dirt traps to simplify cleaning and hygiene management. When securely installed to a suitable, solid wall using appropriate fixings, the entire Prima® Grab Bar range is rated for a maximum user weight of 160 kg.

Intended use

Intended Purpose

To provide a stable, wall-mounted support aid that serves four basic safety functions: pushing up, pulling up, steadying while sitting down, and providing security while transferring position.

Indications for Use

For individuals who may be unsteady, have reduced mobility, or would benefit from extra physical support and leverage to maintain safety and independence in the home.

Intended users

Individuals with reduced mobility, physical weakness, or temporary/permanent disabilities, as well as caregivers assisting them.

Intended Environment

Domestic or residential care environments, suitable for all indoor locations including bathrooms, provided the device is securely fixed to a suitable solid wall.

Contraindication

The device should not be used in the following situations:

- Installation on weak, unsuitable, or unsupported wall surfaces (e.g., standard plasterboard without appropriate structural reinforcement or specialist fixings).
- Use by individuals whose weight exceeds the maximum safe working load of 160 kg.
- Use before verifying that the rail is securely fitted to the wall, or if the end caps are not properly located.
- Use as a step, seat, or for bearing dynamic loads unrelated to human positional transfer and support.

Risk of injury

Please do not use the device and seek help from your local dealer, healthcare professional or provider if any of the following apply:

- The grab bar, concealed fixings, or wall mounts appear damaged, loose, cracked, or unstable.
- The device moves, shifts, or pulls away from the wall when force is applied.
- The ribbed soft grip surface shows signs of significant wear or deterioration that could compromise a secure hold.
- The grip components show signs of deterioration, hardening, splitting, or loss of grip.

Compliance

The manufacturer Gordon Ellis & Co. have a full quality management system accredited by the British Standards Institution to ISO 9001:2015.

EU Declaration of Conformity

This device meets the requirements of the Medical Device Regulation 2017/745, as a Class 1 medical device. The Declaration of Conformity document can be viewed online here: <https://www.gordonellishealthcare.co.uk/prima-grab-bars-2>

Any serious incident arising from use of the device should be reported in a timely manner to the European Authorised Representative (for EU customers) or directly to the manufacturer.

Device testing standards

This device has been tested and meets the requirements of:

- BS EN 12182:2012 - Assistive products for persons with disability. General requirements and test methods
- EN ISO 14971:2019 - Application of risk management to medical devices (risk-assessed)

Expected service life

The expected service life of this product is 5 years, assuming daily use and correct installation, cleaning, and maintenance in accordance with this document.

Storage:

Store indoors in a clean, dry environment ideally between 10°C and 30°C.

Avoid:

- High humidity (above 60%)
- Freezing temperatures
- Excessive heat
- Prolonged direct sunlight
- Damp storage areas

If stored for more than 6 months, the device should be re-inspected before re-issue.

Limitation of liability, and customisation

The manufacturer accepts no liability for damage or injury arising from:

- Customisation or modification: Any alteration to the device - including structural changes, use of non-original spare parts, or unauthorised repairs - renders the device non-compliant and may compromise user safety.
- Improper use: Including use beyond the intended purpose, use on unstable or slippery surfaces, use with unsuitable containers, or applying excessive force.
- Reasonable deterioration of components consistent with expected daily use.

 In the event of contravening point a. or b. the device will no longer retain its CE marking.

Safety information and warnings

Do not use the device until these instructions have been fully read and understood.

If anything is unclear, contact your dealer or healthcare professional before use.

Risk of injury

- Do not use if any part appears damaged, cracked, loose, worn, or unstable.
- Do not exceed the maximum user weight indicated on the label (160 kg)
- Do not use the device as a step, ladder, seating aid, or for bearing dynamic loads unrelated to positional support.
- Do not make any modifications or alterations to the device's structural components.
- Use only original spare parts and accessories supplied or approved by the manufacturer, and do not attempt repairs yourself.

Instructions, device, and label symbols

Instruction symbols: (Image reference: 1a.)

	Important information, useful advice.
	Warning, limitation, or caution.
	Maximum user weight

Label symbols: (Image reference: 1c.)

1	Unique Device Identification (UDI)
2	Manufacturer
3	Country of manufacture
4	Medical device
5	European Conformity mark
6	UK Conformity Assessed mark
7	Lot no. / production date (YY-MM-DD)
8	European Authorised Representative
9	Triman Symbol (France only)
10	European Article Number (EAN)
11	Product linear code-128
12	Product description
13	Product reference number
14	Instructions for use symbol

 Retain the label for future reference.

Device components

(Image reference: 1b.)

Ref.	Component
1	Straight grab rail, available in 12" (300 mm), 16" (400 mm), 18" (450 mm), and 24" (600 mm)
2	Angled grab rail, available in 13" (330 mm) and 16" (400 mm)
3	End cap

Assembly and installation (pg. 3, section 3)

Before assembly, ensure all components are present and that no part is damaged or missing. If any defect is identified, do not use the device and contact your dealer or healthcare provider.

Assembly Instructions

1. Determine orientation:

- Straight rails: Must be fitted either horizontally or vertically.
 - Angled rails: Must be fitted with one section positioned horizontally to allow the user extra leverage.
2. If in doubt, consult a professional tradesperson
 3. Mark fixings: Position the rail against a flat wall and accurately mark the 4 screw positions.
 4. Secure to the wall:
 - The rail may be fitted to most types of solid wall using No. 10 x 50 mm (3/16" x 2") roundhead woodscrews with suitable wall plugs.
 - Depending upon the specific application, other fittings may be suitable; if fitting to plasterboard, ensure suitable specialist fittings are used.
 5. Attach end caps: Once screwed securely into the wall, locate the end caps to the rail and push down to fit them into place.
 6. Check for stability before first use.
 7. See "Cleaning, maintenance, re-issuing" section

Safety notes (Assembly & Installation)

- Wall fixings (screws and wall plugs) are not included; the structural integrity and weight limit of the grab bar are entirely dependent on the competency of the installation and the suitability of the wall structure.
- Always ensure the grab bar is securely fixed to a suitable solid wall. If fitting to plasterboard, you must ensure suitable specialist fittings and structural reinforcements are used.
- If in doubt regarding wall suitability or installation methods, consult a qualified tradesperson.
- Do not install the push-down end caps until the rail has been fully screwed to the wall.
- Before first use, thoroughly check that the rail is securely fitted to the wall.
- The grab bar should be periodically checked for a tight, secure fit - retighten fixings as necessary to ensure continued stability.

Installer / user responsibility

It is the responsibility of the user or installer to:

- Carry out an initial local risk assessment
- Ensure the device is stable before use
- Re-check after cleaning, relocation, or re-issue
- Remove the device from service if instability or wear is detected

Technical Information

Device Materials: (Image reference: 1b.)

1, 2	PP (rail), TPE (grips), Aluminium (core)
3	PP

This device and component parts are manufactured from corrosion-resistant materials, and do not contain natural rubber or latex.

Device dimensions: (Table reference: 2.)

#	Description
1	Clearance from wall to rail handle
2	Rail grip length

#	Description
3	Rail handle width, depth
4	Grab bar total length, width, height
5	Wall fixing hole dimensions, hole size

Cleaning, maintenance, re-issuing, and disposal

To maintain hygiene and preserve the service life of the device, follow the cleaning instructions below.

i DO NOT use abrasive cleaners, wire brushes, scouring pads, or strong chemical disinfectants (e.g., bleach-based or solvent-based cleaners such as acetone or toluene), as these may damage the rail, or grip surfaces.

i ONLY use mild soap or a neutral, non-abrasive disinfectant cleaner (pH 5–9) suitable for plastic and metal surfaces. Always follow the cleaner manufacturer's instructions.

i AVOID prolonged soaking, submersion, or exposure to high-pressure jets of water.

Failure to follow these instructions may lead to corrosion, surface degradation, or reduced safety of the product.

Recommended Cleaning Process

Routine Daily Cleaning

(for general hygiene, household use):

1. Wipe all accessible surfaces of the grab bar with a damp cloth or sponge using mild soap and warm water.
2. Rinse with clean water to remove any detergent residue.
3. Dry thoroughly using a clean, lint-free cloth, paying particular attention to the wall joints and concealed fixing areas.

Enhanced Disinfection

(for enhanced hygiene or after contamination):

1. Use a non-abrasive cold-chemical disinfectant cleaner labeled safe for use on plastics, such as a 70% Isopropyl Alcohol (IPA) solution or a neutral disinfectant spray.
2. Apply the cleaner to the surface and let it sit for the recommended contact time.
3. Wipe clean with a damp cloth and dry fully before use.

Maintenance and Inspection Checklist

Before each use and during periodic inspection, ensure the following:

Inspection Point	Check For	Action if fault found
Grab Bar	Cracks, bends, deformation, or structural damage to the plastic exterior or internal core.	Remove from use immediately and replace.
Soft Grip Surface	Lifting, severe wear, cracking, or loss of grip on the ribbed inserts.	Remove from use immediately and replace.
Wall Fixings & End Caps	Loose screws, wall pull-out, missing end caps, or overall instability when force is applied.	Retighten fixings or consult a tradesperson regarding wall suitability.
Identification Label	Missing or unreadable label.	Contact manufacturer or dealer for replacement.

Routine inspection should be carried out at least every 6 months, or before re-issue. Record the inspection date and findings if required by local procedures.

Re-issuing:

This device is suitable for re-issue, provided

a full inspection and high-level cleaning have been completed in accordance with the above instructions, and the below is met:

i A full inspection has been carried out using the checklist above

i High-level cleaning and disinfection have been completed

i Any worn or damaged components have been replaced using original spare parts

If any component fails inspection and suitable replacement parts are not available, **do not re-issue the device.**

Disposal:

Where possible, disassemble the device into its material groups to support recycling.

This device can be disposed of at many recycling facilities, please check the following for appropriate options in your locality:

- Check your local authority recycling guidelines for plastic goods or mixed materials.
- Separate materials where possible to aid recycling.

i If the device has been used in a healthcare or clinical environment and is considered contaminated, it must be disposed of according to your region's clinical waste regulations.

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