ADJUSTABLE CHAIR RAISER

THE ORIGINAL LANGHAM RAISER, FOR FIRESIDE CHAIRS.



KEY FEATURES

- Quick and easy to fit or remove
- Blocks included to adjust height
- Linked raise for best possible stability
- Recycles well. Long lasting

The original Langham Adjustable Chair Raiser has been tried, tested and trusted by Occupational Therapists for decades. The NEW black plastic replaces the MDF link arms and grey plastic, for a more attractive finish which blends into your household furniture.

The linked design gives the highest level of stability, supplied fully assembled in an easy to fit pack. The load-bearing point is at the bottom of the well, so they are ideal for thinner legs. It is not load bearing on the top.

Height adjusting insert blocks are included to change the raise height from 2-4 inches (5-10cm) depending on the user's needs.

IDEAL FOR: Armchairs.

Do not move the furniture once raisers are fitted.

LP/JLP125PI

Adjustable Chair Raiser



500KG / 78 STONE / 1100 LBS



ADJUSTABLE BED RAISER

DESIGNED FOR RAISING BEDS WITH THIN LEGS.



KFY FFATURES

- Quickly raises beds
- Strong and steady in use
- Versions available for single and double beds
- Recycles well. Long lasting

The Langham Adjustable Linked Bed Raiser can be used to raise beds with legs up to 2 $\frac{1}{2}$ x 3 inch in size (6.5 x 7.5cm). The NEW black plastic replaces the MDF link arms and grey plastic, for a more attractive finish which blends into your household furniture.

It has been carefully designed for beds supported on legs. Each unit is linked across the width of the bed to provide greater strength, in addition, each base is elongated to improve to size of the footprint, giving greater stability and reassurance for the user:

The raise can be altered using the insert blocks under the bed legs to provide three different height adjustments $-2\frac{1}{2}$, $3\frac{1}{2}$ and $4\frac{1}{2}$ inch (6.5, 8.8, 11.5 cm). Legs should only make contact at the bottom of the well. It is not load bearing on top.

IDEAL FOR: Beds with thin legs.

Do not move the furniture once raisers are fitted.

LP/NLP129 Double Beds
LP/NLP129S Single Beds



500KG / 78 STONE / 1100 LBS